


# Imaging system manual

Written by Itamar Shachen Tov & Elad Vizenblit

## Requirements

For a new computer – in the sound\_server folder, run in the cmd: “npm i” this will download all the necessary libraries for running the program. It will save them in a folder called node\_modules

```
C:\Users\user\Desktop\imaging_system\ViRMEn 2016-02-12\sound_server>npm i
```




Create a folder called “log” inside the virmen folder.

(C:\Users\user\Desktop\imaging\_system\ViRMEn 2016-02-12\log)

Create a folder called “config” inside the virmen folder.


(C:\Users\user\Desktop\imaging\_system\ViRMEn 2016-02-12\config)



## For running the program:

Open Matlab 2018b



Go to the “C:\Users\user\Desktop\imaging\_system\ViRMEn 2016-02-12” folder in the matlab file explorer.




This is the way you have to write virmen in the terminal followed by Enter



You will get this window





### What do we see in the image above?


On the right side this we have 4 world implemented. (The experiment use only 3 of them as we know- chess, stripes, black room). Please don't touch it.

MoveWithDAQ is our movement function. This is the one we have implemented for this experiment. The others will not work correctly.

The bar at the top of the screen is used for creating a new worlds. Please don't click there as well. There is no reason. (The only button necessary is the run button, it is being used in the next step)


Back to our how to run manual:

after setting everything, run the session



you will enter the configuration of the session.

After pressing "select" It will be saved in the log folder, in a file called "config.json" The default configuration file when getting to this screen is the last configuration used.



- 1) Number of trials for the session
- 2) Percentage of time running in the goal range in small reward world in order to open the valve
- 3) Reward size- how much time to open the valve for the small reward.
- 4) Which type of sound policy to activate.
- 5) Choose the hallway that will represent the large reward hallway.
- 6) How much time (in sec) to be in the reward room.
- 7) How much time (in sec) to be in the black room.
- 8) Target speed (cm/s)
- 9) How much deviation from the target speed we allowed to each side(cm/s). (i.e it sets to 20-40 in this configuration)
- 10) in out of range scenario, how much deviation between each step
- 11) As of 2 – but here it is for large reward room.
- 12) As of 3 – but here it is for large reard room
- 13) The length of the trace (in cm, not precise)

## Using a previously used configuration setting


For using a config.json file from other session, copy its content and paste it in the default config file in config\config.json (you can use CTRL+C CTRL+V).

- If by any reason there is no config.json in config folder. It will create a default file when starting a session.

## results

after running the session, all the data will be saved in the log folder, under the session's date and time

There are 5 files in the folder:



## How the data is structured

A-B\_Leakport\_record (4Xn)

timestamp	...
A state(0-4 Volt)	...
B state(0-4 Volt)	...
Lickport state	...

Reward(2Xn)

Timestamps	...
Reward gaved	...

Velocity (3Xn)

timestamp	...
Velocity (cm\s)	...
Current leakport status	...

Synchronization signal (2Xn)

timestamp	...
Random signals	...


Pay attention - the Synchronization file contains timestamps for the first 200 samples each writing (when we start the trial and every time we get to the end of the trace). The rest are filled with zeros.

### **How to set sound files:**

The sound files are in this file path:

C:\Users\user\Desktop\imaging\_system\ViRMEn 2016-02-12\sound\_server

And these are the sound files the program use.




sound1000.wav  
sound1200.wav  
sound1400.wav  
sound1600.wav  
sound1800.wav  
sound2000.wav  
sound2200.wav  
sound2400.wav  
sound2600.wav  
sound2800.wav  
sound3000.wav  
sound3200.wav  
sound3400.wav  
sound3600.wav  
sound3800.wav  
sound4000.wav  
sound4200.wav  
sound4400.wav  
sound4600.wav  
sound4800.wav  
sound5000.wav

The desired sound is sound3000.wav for all auditory feedbacks.

### Preferred speed:

Sound3000.wav will be played when it runs in desired velocity range


In case it doesn't run in desired velocity range- no sound file will be played.



### Gradual symmetrical tones:

Sound3000.wav will be played when it runs in desired velocity range

In case it goes out of range – every step we go out from range in both sides the sound file number will be increased by 200.




### Gradual asymmetrical tones:

Sound3000.wav will be played when it runs in desired velocity range

In case it goes out of range, it splits into two – too slow and too fast.

For every step we go too fast, the sound file number is increased by 200.

For every step we go too slow, the sound file number is decreased by 200.





So all you have to do is replace the sound files with the desired sounds you want to have based on this configuration.